

InSight®

qLabs® Veterinary Coagulation Analyser

Calibrated and validated specifically
for use in dogs and cats



COAGULATION SCREEN. ONE DROP WHOLE BLOOD.

Tel: 0300 303 3509 Email: info.vpg@synlab.co.uk Web: vet.synlab.co.uk

VPG Exeter, Unit 1B Exeter Science Park,
6 Babbage Way, Clyst Honiton, Exeter, EX5 2FN



qLabs® Veterinary Coagulation Analyser

In-Practice Point of Care Coagulation Testing

Calibrated and validated specifically for use in dogs and cats

Key Features

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/aPTT tested simultaneously on 1 test strip
- Room temperature test strip storage
- Fill test strip directly from syringe
- QC material available
- No need to transfer sample into citrate tube
- 'Test count up' for immediate assessment

Simple Test Procedure

1. Insert test strip
2. Apply a drop of blood
3. Results displayed on LCD screen

Helping Diagnose Conditions in Relation to:

- Pre-operative checks for at risk animals
- Toxins and systemic diseases
- Hereditary bleeding disorders
- Excessive bleeding
- Disseminated Intravascular Coagulopathy (DIC)

Cartridge Types

- PT/aPTT combo cartridge (12)

Kit Contents

- qLabs® Coagulation Analyser
- Power Adapter
- Lancets (10)
- Lancing Device
- User Manual
- Carry Case



2 Year Warranty